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Tel +39 0185 939481 - Fax +39 0185 938595

info@larident.it - www.larident.it

# MANUFACTURER'S INSTRUCTIONS FOR THE TREATMENT OF MEDICAL DEVICES ACCORDING TO THE EN ISO 17664:2017 STANDARDS

**MANUFACTURER:** Larident Srl

**REF:** C90/C92/C99/C95

---

**WARNING NOTICE:** THE PRODUCT IS NOT SUPPLIED STERILE.  
STERILIZE IT BEFORE USE

---

**Any serious incident related to this medical device must be reported to the manufacturer and the competent authorities.**

**Restriction on treatments:** the product is for single use

**Area of use:** support and cleaning of endodontic instruments

**Storage and transport:** the product does not have specific needs

**Method of use:** clean the file or reamer files by rotating the instrument, guided by a hand inside the sponge. The impurities left in the endodontic instruments will be removed by the porosity of the sponges

**Disinfection and sterilization:** sterilize in an autoclave machine before use at 134° for 4 minutes. The product could slightly change colour but this will not change its usage or performance

**Maintenance:** no specific requirement

**Disposal:** dispose of the product according to the correct disposal procedures

**Function check and test:** visual check for any damage or deformation of the product

**Storage:** the product does not have specific needs

The instructions listed above have been validated by the manufacturer of the medical devices as suitable for preparing the device for use.

It is the user's responsibility to ensure that the preparation process, the equipment the materials used and the staff obtain the desired results.

Any deviation from the recommended instructions must be carefully checked with regards to its efficiency and possible adverse consequences.

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**LARIDENT S.r.l.**