VOCO Clip Flow

CE 0482

EN Instructions for use
MD EU Medical device

Product description:

Clip Flow is a flowable, light-curing material for temporary fillings, sealings and treatments. Thanks to its elastic consistency, the material is easily removable, also in case of undercuts. Therefore, post-treatment of the cavity is not necessary.

Indications:

- Temporary inlay and onlay treatments of the cavity
- Sealing of openings for implant screws
- Relining material for temporary crowns and bridges
- Block-out material for retentive areas in the dental arch, e.g. before taking impressions
- Covering of the gingival margin
- Fixing of resin matrix during filling placement
- Temporary filling of cavities and sealing of root canal orifices

Contraindications:

Clip Flow contains methacrylate and BHT. **Clip Flow** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Clip Flow is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Clip Flow should only be used by a professionally trained dental practitioner.

Use:

Preparation:

If necessary, carry out preparation, lining and impression. Rinse and dry the cavity before application of **Clip Flow**.

Application:

Apply **Clip Flow** directly from the NDT® syringe using the respective cannula (type 46). In order to ensure the functioning of the non-dripping syringe do not pull back the plunger during or after use.

If necessary, spread **Clip Flow** with an adequate instrument (e.g. probe). Light-cure with a polymerisation light (minimum 500 mW/cm²) a film thickness up to 1 mm for 10 s, (e.g. in case of covering the gingival margin), a thickness up to 5 mm for 20 s and more than 5 mm for 40 s. Hold the light emission tip of the device as close as possible to the surface of the filling. Otherwise, the curing depth may be reduced.

Finishing:

Excess material can be removed with a scalpel or round bur. Silicone points are suitable for finishing.

Removal of the temporary filling:

Insert an appropriate instrument (e.g. probe or curette) sideways into the temporary filling and extract **Clip Flow** from the cavity.

Warnings, precautionary measures:

- Adhesion to dental resins or linings based on (meth)acrylate may occur. In such cases apply a suitable insulating agent (e.g. glycerine, vaseline).
- In case of eye contact, rinse eyes immediately with plenty of water and seek medical advice.
- The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Clip Flow**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Clip Flow**.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Clip Flow** must be reported to VOCO GmbH and the responsible authority.

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