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# AFFINIS<sup>®</sup> putty soft / putty super soft / fast putty soft

## Definition

AFFINIS is a silicone-based impression material for use in dentistry.

## Material type

Polyvinylsiloxane, addition-type silicone elastomer.

### putty soft

Base: brown  
Catalyst: light grey  
– ISO 4823, Type 0, putty consistency

### putty super soft

Base: orange  
Catalyst: light grey  
– ISO 4823, Type 0, putty consistency

### fast putty soft

Base: brown  
Catalyst: light grey  
– ISO 4823, Type 0, putty consistency


## Indications

putty soft / putty super soft / fast putty soft

- Primary impression in the putty-wash impression technique
- Tray material for the simultaneous mixing technique
- Tray material for the copperband impression technique
- Impressions for study models, orthodontic models, matrices

## Clinical time

### AFFINIS putty soft

### AFFINIS putty super soft

### AFFINIS fast putty soft

## Important

- Check intraorally if the material is completely set before removing.
- High temperatures will speed up the setting process and low temperatures will slow it down.
- For extended working time, material can be placed in the refrigerator prior to use.

## Contra-indications

Provided that the product is used as directed, there are no known contra-indications.

## Side effects / Interactions

Polyvinylsiloxanes have very good biological compatibility and up to now no harmful reactions or secondary effects on patients and/or dental personnel are known. Impression materials are intended to set in the mouth of the patient, however, they should remain in the mouth not more than twice the setting time. Although they have reasonably high tear strength, care should be taken that no portion of the impression material remains in the interdental spaces or in the sulcus. Undercuts should, in certain instances, be blocked out before taking the impression.

## Trays

Tray selection depends on impression technique and professional preference of stock trays, double arch bite trays or custom trays. For perfect adhesion, we recommend applying a thin layer of Coltène<sup>®</sup> Adhesive or any other brand of adhesive specified for use with polyvinylsiloxane impression materials.

## Dosage

Measure out level scoopfulls of Base and Catalyst. **Attention:** Avoid interchanging the color coded measuring-scoops, each belonging to one of the components, since there is a risk of contamination!

## Mixing

1. Take the base material and the catalyst with the fingertips and start mixing. (Always wear gloves)
2. Knead the material with the fingertips until the color is uniform.

**Attention:** Only use the Base and Catalyst from the same package.

## Recommendations for use

### «double mix technique»

While mixing, the dentist may begin to apply the wash-material around the preparation. Tray loading and application of wash-material around the preparation have to be completed at the same time. Keep the oral tip immersed in the material at all times to avoid trapping air. *Immediately place tray filled with heavy body material into the patient's mouth. Press for 2–3 s and hold in position until the material is completely set.*

### «two step technique»

If a 2-Step «putty/wash» technique is used, the primary impression made with putty material, must be carefully cleaned and dried before use. This will help guarantee good adhesion between the wash (low viscosity material) and putty tray material.

**For extended working time, material can be placed in the refrigerator prior to use.**

## Important

Always wear gloves. Some latex gloves and surfaces contaminated by these gloves (teeth, preparations, retraction cords, etc.) may interfere with the setting process of polyvinylsiloxanes. The material and surfaces where the impression will take place (teeth, preparations, retraction cords, etc.) should only come into contact with thoroughly washed and rinsed gloves (wash 15 s with a detergent, rinse with lukewarm tap water for another 15 s) or vinyl-/nitrile gloves. To determine if your gloves will interfere with the setting process, extrude a small amount of material onto the palm of a gloved hand. Close your hand around the material and see if it sets within the time limits listed above. Compounds containing eugenol or hemostatic substances can also impede perfect setting. If (H<sub>2</sub>O<sub>2</sub>) hydrogen peroxide is used for disinfection, it is recommended to thoroughly rinse with lukewarm water afterwards in order to avoid bubble formation. Intensive mixing or high temperatures speed up, low temperatures slow down the setting.

**Check intraorally if the material is completely set before removing from the mouth.**

## Disinfection

The impression should be rinsed under (cold) running tap water after removal from the mouth. After rinsing, disinfection with a suitable commercial dental disinfectant solution will not affect the impression surface or dimensions. Acrylic trays should be protected against water absorption.

## Optional

The following points should be observed when auto-

## claving impressions:

1. Only use autoclavable components (e.g. PRESIDENT Tray AC, Coltène Adhesive AC).
2. Rinse and clean impression thoroughly under flowing, lukewarm water.
3. Impressions can be autoclaved directly after taking the impression.
4. Autoclave in steam steriliser at 134°C/ 273°F (Prior program)

## Important

When sterilising implant impressions, be sure clarify in advance with the manufacturer, whether or not the implant components (e.g. impression posts, etc.) are autoclavable.

## Manufacture of models

Die stone models may be poured after 30 min at the earliest. The impression remains dimensionally stable for a practically unlimited period of time (at least 7 days). The surface tension will be reduced and pouring will be facilitated if the impression is briefly washed out with a detergent and rinsed thoroughly in luke-warm clear water afterwards. Coltène<sup>®</sup> impression materials are compatible with standardized dental stones (e.g. Fuji-rock, Vel-mix), epoxy resin and polyurethane.

## Electro plating

Compatible with commercially available copper plating and silver plating baths.

## Cleaning of trays

The impression can be removed with a blunt instrument. Soaking in a universal commercial solvent or light petrol will dissolve the Coltène<sup>®</sup> Adhesive. These solvents should only be used in a well-ventilated area. The trays can be cleaned and disinfected as usual.

## Shelf life and storage

AFFINIS fulfils the intended purpose at least until the expiry date, which is displayed on the containers, in well sealed containers at temperatures 15–23 °C / 59–73 °F and 50 % relative air humidity. Impressions should be stored at normal room temperature, avoid exposure to heat and sun.

## Marking

The expiry date and number [LOT] are shown on the package and cartridge.

## Technical data ISO 4823:2000

Measurements are based on 23 °C / 73 °F room temperature and 50 % relative humidity.

### putty soft

Mixing time (15ml): 0:30 min  
Time in mouth (35 °C / 95 °F): 2:00 min

### putty super soft

Mixing time (15ml): 0:25 min  
Time in mouth (35 °C / 95 °F): 2:00 min

### fast putty soft

Mixing time (15ml): 0:30 min  
Time in mouth (35 °C / 95 °F): 1:50 min

## Date of issue

03-2014

## Caution:

Federal law restricts this device to sale by or on the order of a dentist.

Only supplied to dentists and dental laboratories or upon their instructions.

For MSDS see: [www.coltene.com](http://www.coltene.com)

# AFFINIS<sup>®</sup> putty soft / putty super soft / fast putty soft

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**Coltène/Whaledent AG** 🏢