Manufactured by GDF GmbH, Rosbach | Germany
for Brilliant Smile Sweden AB | Backa Bergögata 16 | 422 46 Hisings Backa | Sweden

Min. 3° C - Max. 9° C

CE0297

Refrigerate the bleaching gel at 3-9°C. Do not freeze.

Product contains hydrogen peroxide, potassium nitrate, sodium fluoride, sodium hydroxide, glycol derivatives, thickening agent and dye.

(EN)

Treatment periods

PowerWhite^{pro} 3 x 10 min

Directions for use

PowerWhite^{pro} is a bleaching gel intended for medical in-office use only. The material is used by the dentist to brighten individual teeth. The gel is ideally suited for brightening devital teeth including endodontically treated teeth. The 4:1 auto-mix syringe (cartridge) contains bleaching gel and activator. One syringe contains 5 g of material. The hydrogen peroxide content of the ready-mixed gel is 35%.

Indications

- Brightening of discoloured, devital teeth
- Brightening of discoloured teeth caused by disease, by injury or iatrogenically, as an alternative to a crown or veneer
- Treatment of tooth discolourations caused by prescription drugs (e. g. tetracycline)

Contraindications

- Exposed dentine
- Defective filling margins
- Hypersensitive patients
- Allergy to or intolerance of the ingredients
- Pregnancy and lactation (due to insufficient data on the effects on pregnant and lactating women)

Many patients may experience slight sensitisation in their teeth during or shortly after treatment, but this normally subsides after a short time.

Stop treatment immediately if undue pain is being experienced.

Preparation

Determine the shade of the tooth by comparison with a Vita© shade guide (photo 1). For assessment, re-arrange the Vita© shade guide as listed below.

Values arranged in accordance with the shades:

A1 A2 A3 A3.5 B1 B2 B3 B4 C1 C2 C3 C4 D1 D2 D3 D4

Values arranged in accordance with the brightness:

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B2 to B1 = increase of value by 2

A3 to A2 = increase of value by 4

Insert a cheek or lip retractor (photo 2). In case you start by cleaning the patient's teeth with air abrasion, make sure the powder is bicarbonate based. Polish teeth thoroughly with pumice without additives (photo 3). Rinse and blow completely dry. Do not let the patient rinse their mouth during treatment, as the glucoprotein in saliva significantly decreases the effect of the bleaching gel. Apply the gingiva protector along the gum line (photo 4), and cure it with the BrilliantSmile LED Plasma Light (photo 5). Do not anaesthetise the patient.

The ready-mixed bleaching gel irritates the skin and the mucous membranes, and it causes severe eye damage. Always wear protective gloves and goggles when using the material! Avoid contact of the material with eyes, skin and mucous membranes.

Application

Remove the cap from the cartridge and replace with the mixing tip. Press the gel out of the cartridge until a uniformly colored gel emerges from the mixing tip. Discard the inhomogeneous material which emerged initially.

Storage at higher temperatures for a longer period of time can cause the development of isolated bubbles within the material. This may lead to inhomogeneous, liquid mixtures, especially when there is less than a quarter of the material left inside the double cartridge. In this case, do not use the remaining material.

Apply the gel in a layer of 1-2 mm thickness onto the labial surfaces of the teeth to be bleached (photo 6). The mixing tip is very well suited for spreading the gel. Any gel on the gums or soft tissue must be removed immediately. The effectiveness of the product is increased by use of the BrilliantSmile LED Plasma Light (photo 7).

The patient and dentist should wear safety goggles.

Remove the gel by suction after 10 min, rinse off any residue with water and blow dry (photo 8). The process is repeated two more times. Before repeating the process, make sure the gingiva protection remains intact.

After the last bleaching process, remove the gel by suction, and thoroughly rinse off any residue with water and blow dry. Determine the tooth shade once again (photo 10). Afterward, the teeth should be remineralized with a fluoride solution or gel. Such solution or gel is also suitable for treatment of sensitization or pain occurring during the bleaching process. After treatment, the teeth are particularly susceptible to discoloration. Offer the patient the BrilliantSmile NanoSeal Total+ varnish (photo 9), otherwise the patient must avoid acidic and discoloring foods, beverages and tobacco for 72 hours, in addition to rinsing with fluoride twice daily, in order to keep the best possible result.

If the material remaining in the cartridge is sufficient for another treatment, remove the mixing tip and close the cartridge with the blue cap. Store the material in a refrigerator. Use a new mixing tip for each patient.

Note: Composite fillings or ceramic and metal restorations cannot be bleached. Such fillings or restorations should be replaced at least two weeks after the treatment.

Devital or endodontically treated teeth

The gel is also applied through an access opening on the lingual side of the tooth. First, check the height of the alveolar ridge. Then completely remove the filling material from the pulp chamber.

Only apply one layer of glass ionomer filling with a thickness of 1- 2 mm to the bottom of the pulp chamber on the same level with the alveolar ridge or coronally. First apply the gingiva protection to both the labial and the lingual sides of the tooth. Now, apply the bleaching gel to the labial surface as described under "Application" above. Additionally, apply the gel through the access opening. Let the gel act for 10 min. Afterwards, carefully aspirate the gel and lightly rinse with water. If the bleaching effect is insufficient, this procedure can also be repeated twice until the desired effect has been achieved. After the final bleaching procedure, aspirate the gel and rinse thoroughly with water. Carefully remove the gingiva protection. Next, close the access opening with a temporary filling material which does not contain eugenol. Wait two weeks before preparing the permanent filling.

Shelf-life/Storage

Gel should not be used after its expiration date, as a sufficient bleaching performance can then no longer be guaranteed.

Possible sources of error

Insufficient bleaching effect	Repeat bleaching process (up to four treatment cycles in one sitting). Additional bleaching treatments can be performed after one week.
Patient is experiencing pain	Gel is in contact with exposed dentine or a loose filling (cover before treatment) Patient has particularly sensitive teeth (depending on the individual tooth structure)
White spots on the enamel after bleaching	Local demineralization or dehydration. Reversible. Disappears after fluoridation.
Gingiva is coloured white after bleaching	The gingiva is not adequately protected, or the protection is loose.
The gel from the cartridge has inhomogeneous colour (colourless/red)	The gel has air bubbles in it as a result of being stored at too high a temperature or for too long.
Reduced composite-enamel adhesion after bleaching procedure	Peroxide restoration. Leave for 1-2 weeks before restoration.

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