Prior to use, carefully read the instructions for use. **EN**

Miracle Mix™ CAPSULE

RADIOPAQUE SILVER ALLOY-GLASS IONOMER CEMENT RESTORATIVE IN CAPSULES

For use only by a dental professional in the recommended indications.

INDICATIONS FOR USE

1. Core build-up
2. Restoration of primary teeth
3. Cases where radiopacity is required
4. As a base or liner
5. Class I, limited Class II, temporary fillings, etc.

CONTRAINDICATIONS

1. Direct pulp capping.
2. Avoid use of this product in patients with known allergies to glass ionomer cement, silver, tin or copper.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	4.8 / 1.0
Mixing Time (sec.)	10"
Time of initial extrusion after mixing	10" maximum
Working Time (23°C/73°F) (min., sec.)	1'00"
Net Setting Time (37°C/99°F) (min., sec.)	4'50"
Finishing Time (37°C/99°F, from start of mixing) (min., sec.)	5'00"

ISO 9917-1 Glass polyalkenoate cement, Restoration

1. CAVITY PREPARATION

1. Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.

- b) Wash and dry but DO NOT DESICCATE. CAVITY CONDITIONER (10 seconds) or DENTIN CONDITIONER (20 seconds) is recommended to remove the smear layer (Fig. 1).

2. CAPSULE ACTIVATION / MIXING

- a) Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
- b) To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds (Fig. 3).

Note:

- 1) Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
- 2) The capsule should be activated just before mixing and used immediately.

- 3) Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000 RPM) (Fig. 4).

3. RESTORATIVE TECHNIQUE

- a) Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER. Time of initial extrusion after mixing: 10 seconds maximum.

- b) Click twice to prime the capsule, then dispense. (Fig. 5). The working time is 1 minute from start of mixing at 23°C (73°F). Higher temperature will shorten working time.

- c) Remove surface moisture but DO NOT DESICCATE.
- d) Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation (Fig. 6). Avoid air bubbles. Note:

- 1) To adjust the direction of the nozzle, hold the applicator with the capsule towards you and twist.
- 2) To remove the used capsule, push the applicer leverage button. Twist the capsule and pull upwards.

- e) Form the contour.
- f) Setting time is 5'00" from start of mixing.

4. FINISHING & POLISHING
Use standard techniques (Fig. 7).

STORAGE

Recommended for optimal performance, store at temperature of 4-25°C (39.2-77.0°F).

PACKAGE

Package: 50 capsules
Average contents per capsule: 0.52 g silver alloy ionomer cement and 0.11 g (0.09 mL) liquid
Minimum net volume of mixed cement per capsule: 0.09 mL
Option: GC CAPSULE APPLIER (1 piece)

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
4. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com>

They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS
United Kingdom

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