

Implant Compatibility List for NeoBiotech Implant Removal Kit

NOTE: For other Manufacturers not listed, the clinician is advised to complete a "try in" by selecting a corresponding analog or implant to thread into whereby sequentially starting with the widest diameter FRS screw M2.5 followed by, in descending order, the next available diameter. This will ensure compatibility.

| Corporation | Implant System | Implant Ø Diameter | Screw Size to use | Solution Kit Compatibility | Product Name |
|----------------|--------------------------------------|-----------------------|----------------------|-------------------------------|-----------------|
| Nobel Biocare® | Branemark System® MK III Groovy NP | 3.3mm | M1.6 | YES | FRS 16 |
| | Branemark System MK III Groovy RP | 3.75mm | M2.0 | YES | FRS 20 |
| | Branemark System MK III Groovy RP | 4.0mm | M2.0 | YES | FRS 20 |
| | Branemark System MK III Groovy WP | 5.0mm | M2.5 | YES | FRS 25 |
| | NobelSpeedy® System MK III Groovy NP | 3.3mm | M1.6 | YES | FRS 16 |
| | NobelSpeedy System MK III Groovy RP | 4.0mm | M2.0 | YES | FRS 20 |
| | NobelSpeedy System MK III Groovy WP | 5.0mm | M2.0 | YES | FRS 20 |
| | NobelSpeedy System MK III Groovy WP | 6.0mm | M2.5 | YES | FRS 25 |
| | NobelReplace® Straight Groovy NP | 3.5mm | M1.8 | YES | FRS 18 |
| | NobelReplace Straight Groovy RP | 4.0mm | M2.0 | YES | FRS 20 |
| | NobelReplace Straight Groovy WP | 5.0mm | M2.0 | YES | FRS 20 |
| | NobelSpeedy Replace NP | 3.5mm | M1.8 | YES | FRS 18 |
| | NobelSpeedy Replace RP | 4.0mm | M2.0 | YES | FRS 20 |
| | NobelSpeedy Replace WP | 5.0mm | M2.0 | YES | FRS 20 |
| | NobelSpeedy Replace 6.0 | 6.0mm | M2.0 | YES | FRS 20 |
| | NobelReplace Tapered NP | 3.5mm | M1.8 | YES | FRS 18 |
| | NobelReplace Tapered RP | 4.3mm | M2.0 | YES | FRS 20 |
| | NobelReplace Tapered WP | 5.0mm | M2.0 | YES | FRS 20 |
| | NobelReplace Tapered | 6.0mm | M2.0 | YES | FRS 20 |
| | NobelActive® | 3.5mm | M1.54 | SPECIAL ORDER | FRS 154 |
| Straumann® | Standard | 3.3 RN | M2.0 | YES | FRS 20 |
| | Standard | 4.1 RN | M2.0 | YES | FRS 20 |
| | Standard | 4.8 RN | M2.0 | YES | FRS 20 |
| | Standard | 4.8 WN | M2.0 | YES | FRS 20 |
| | Standard Plus | 3.3 NN | M1.6 or M1.8 | Verify | FRS 16 / FRS 18 |
| | Standard Plus | 3.3 RN | M2.0 | YES | FRS 20 |
| | Standard Plus | 4.1 RN | M2.0 | YES | FRS 20 |
| | Standard Plus | 4.8 RN | M2.0 | YES | FRS 20 |
| | Standard Plus | 4.8 WN | M2.0 | YES | FRS 20 |
| | Tapered Effect | 3.3 RN | M2.0 | YES | FRS 20 |
| | Tapered Effect | 4.1 RN | M2.0 | YES | FRS 20 |
| | Tapered Effect | 4.8 RN | M2.0 | YES | FRS 20 |
| | Bone Level Implant | 3.3 NC | M1.6 or M1.8 | Verify | FRS 16 / FRS 18 |
| | Bone Level Implant | 4.1 NC | M1.6 or M1.8 | Verify | FRS 16 / FRS 18 |
| | Bone Level Implant | 4.8 NC | M1.6 or M1.8 | Verify | FRS 16 / FRS 18 |
| BIOMET 3i | Certain® Internal | 3.4 | M1.6 | YES | FRS 16 |
| | Certain Internal | 4.1 | M1.6 | YES | FRS 16 |
| | Certain Internal | 5.0 | M1.6 | YES | FRS 16 |

Implant Compatibility List for NeoBiotech Implant Removal Kit

NOTE: For other Manufacturers not listed, the clinician is advised to complete a "try in" by selecting a corresponding analog or implant to thread into whereby sequentially starting with the widest diameter FRS screw M2.5 followed by, in descending order, the next available diameter. This will ensure compatibility.

| Corporation | Implant System | Implant Ø Diameter | Screw Size to use | Solution Kit Compatibility | Product Name |
|-------------------------------|---|-----------------------|----------------------|-------------------------------|-----------------|
| | Certain® Internal | 6.0 | M1.6 | YES | FRS 16 |
| | External | 3.4 | M2.0 | YES | FRS 20 |
| | External | 4.1 | M2.0 | YES | FRS 20 |
| | External | 5.0 | M2.0 | YES | FRS 20 |
| | External | 6.0 | M2.0 | YES | FRS 20 |
| Zimmer® (Screw-Vent®) | Screw Advent Implant | 3.3 | 1-72UNF | YES | FRS 172 |
| | Screw Advent Implant | 3.7 | 1-72UNF | YES | FRS 172 |
| | Screw Advent Implant | 4.7 | 1-72UNF | YES | FRS 172 |
| | Tapered Screw Advent Implant | 3.7 | 1-72UNF | YES | FRS 172 |
| | Tapered Screw Advent Implant | 4.7 | 1-72UNF | YES | FRS 172 |
| | Tapered Screw Advent Implant | 6.0 | 1-72UNF | YES | FRS 172 |
| | Advent® Implant | 3.7 | 1-72UNF | YES | FRS 172 |
| | Advent Implant | 4.7 | 1-72UNF | YES | FRS 172 |
| | Advent Implant | 4.7 | 1-72UNF | YES | FRS 172 |
| Zimmer® (Swissplus®) | Straight Swissplus Implant | 4.1 | M2.0 | YES | FRS 20 |
| | Straight Swissplus Implant | 4.8 | M2.0 | YES | FRS 20 |
| | Tapered Swissplus Implant | 3.7 | 1-72UNF | YES | FRS 172 |
| | Tapered Swissplus Implant | 3.7 | M2.0 | YES | FRS 20 |
| | Tapered Swissplus Implant | 4.8 | M2.0 | YES | FRS 20 |
| Zimmer® (Spline®) | Spline Twist Implants | 3.75mmD | 1-72UNF | YES | FRS 172 |
| | Spline Twist Implants | 5.0mmD | 1-72UNF | YES | FRS 172 |
| | Spline Reliance Cylinder Implants | 3.25mmD | 1-72UNF | Verify | FRS 172 |
| | Spline Reliance Cylinder Implants | 4.0mmD | 1-72UNF | YES | FRS 172 |
| | Spline Reliance Cylinder Implants | 5.0mmD | 1-72UNF | YES | FRS 172 |
| Astra Tech® | OsseoSpeed | 3.0S | M1.4 | NO | FRS 14 |
| | OsseoSpeed | 3.5S | M1.6 | YES | FRS 16 |
| | OsseoSpeed | 4.0S | M1.6 | YES | FRS 16 |
| | OsseoSpeed | 4.5 | M2.0 | YES | FRS 20 |
| | OsseoSpeed | 5.0 | M2.0 | YES | FRS 20 |
| | OsseoSpeed | 5.0S | M2.0 | YES | FRS 20 |
| Keystone Dental (Lifecore) | PrimaConnex® Internal Connection Tapered Implants | 3.5 | M1.8 | YES | FRS 18 |
| | PrimaConnex Internal Connection Tapered Implants | 4.1 | M1.8 | YES | FRS 18 |
| | PrimaConnex Internal Connection Tapered Implants | 5.0 | M1.8 | YES | FRS 18 |
| | PrimaConnex Internal Connection Straight Implants | 3.3 | M1.8 | YES | FRS 18 |
| | PrimaConnex Internal Connection Straight Implants | 4.0 | M1.8 | YES | FRS 18 |
| | PrimaConnex Internal Connection Straight Implants | 5.0 | M1.8 | YES | FRS 18 |
| | RENOVA® RBM Tapered Implants | 3.75 | M1.8 | YES | FRS 18 |
| | RENOVA RBM Tapered Implants | 4.5/4.75 | M1.8 | YES | FRS 18 |

Implant Compatibility List for NeoBiotech Implant Removal Kit

NOTE: For other Manufacturers not listed, the clinician is advised to complete a "try in" by selecting a corresponding analog or implant to thread into whereby sequentially starting with the widest diameter FRS screw M2.5 followed by, in descending order, the next available diameter. This will ensure compatibility.

| Corporation | Implant System | Implant Ø Diameter | Screw Size to use | Solution Kit Compatibility | Product Name |
|-------------|--------------------------------------|-----------------------|----------------------|-------------------------------|-----------------|
| | RENOVA RBM Straight Implants | 3.75 | M1.8 | YES | FRS 18 |
| | RENOVA RBM Straight Implants | 4.5/4.75 | M1.8 | YES | FRS 18 |
| | RESTORE® Small Diameter | 3.3 | 1-72UNF | YES | FRS 172 |
| | RESTORE Regular Diameter | 3.75 & 4.0 | M2.0 | YES | FRS 20 |
| | RESTORE Small Diameter Wide Diameter | 5.0 & 6.0 | M2.5 | YES | FRS 25 |
| | STAGE-1® system | 3.3 (RDS) | M2.0 | YES | FRS 20 |
| | STAGE-1 system | 4.1 (RDS) | M2.0 | YES | FRS 20 |
| | STAGE-1 system | 4.8 (RDS) | M2.0 | YES | FRS 20 |
| | STAGE-1 system | 4.8 (WDS) | M2.0 | YES | FRS 20 |
| | STAGE-1 system | 5.5 (WDS) | M2.0 | YES | FRS 20 |
| | STAGE-1 system | 6.3 (WDS) | M2.0 | YES | FRS 20 |
| MIS® | MISTRAL Internal Octagon System | 3.75 | M2.0 | YES | FRS 20 |
| | MISTRAL Internal Octagon System | 4.1 | M2.0 | YES | FRS 20 |
| | MISTRAL Internal Octagon System | 4.8 | M2.0 | YES | FRS 20 |
| | Biocom Internal | 3.3 | 1-72UNF | YES | FRS 172 |
| | Biocom Internal | 3.75 | 1-72UNF | YES | FRS 172 |
| | Biocom Internal | 4.2 | 1-72UNF | YES | FRS 172 |
| | Biocom Internal | 5.0 | 1-72UNF | YES | FRS 172 |
| | Biocom Internal | 6.0 | 1-72UNF | YES | FRS 172 |
| | Seven® Internal | 3.75 | 1-72UNF | YES | FRS 172 |
| | Seven Internal | 4.2 | 1-72UNF | YES | FRS 172 |
| | Seven Internal | 5.0 | 1-72UNF | YES | FRS 172 |
| | Seven Internal | 6.0 | 1-72UNF | YES | FRS 172 |
| | Lance Internal | 3.75 | 1-72UNF | YES | FRS 172 |
| | Lance Internal | 4.2 | 1-72UNF | YES | FRS 172 |
| | Lance Internal | 5.0 | 1-72UNF | YES | FRS 172 |
| Dentsply™ | Ankylos® internal | 3.5 | M1.8 | YES | FRS 18 |
| | Ankylos internal | 4.5 | M1.8 | YES | FRS 18 |
| | Ankylos internal | 5.5 | M1.8 | YES | FRS 18 |
| | Ankylos internal | 7.0 | M1.8 | YES | FRS 18 |
| | Xive® S internal Implant | 3.0 | M1.6 | YES | FRS 16 |
| | Xive S internal Implant | 3.4 | M1.6 | YES | FRS 16 |
| | Xive S internal Implant | 3.8 | M1.6 | YES | FRS 16 |
| | Xive S internal Implant | 4.5 | M1.6 | YES | FRS 16 |
| | Xive S internal Implant | 5.5 | M1.6 | YES | FRS 16 |
| | Xive TG External Implant | 3.4 | M1.4 | Verify | FRS 14 |
| | Xive TG External Implant | 3.8 | M1.4 | Verify | FRS 14 |
| | Xive TG External Implant | 4.5 | M1.4 | Verify | FRS 14 |

Implant Compatibility List for NeoBiotech Implant Removal Kit

NOTE: For other Manufacturers not listed, the clinician is advised to complete a "try in" by selecting a corresponding analog or implant to thread into whereby sequentially starting with the widest diameter FRS screw M2.5 followed by, in descending order, the next available diameter. This will ensure compatibility.

| Corporation | Implant System | Implant Ø Diameter | Screw Size to use | Solution Kit Compatibility | Product Name |
|--------------|----------------------------|--------------------|-------------------|----------------------------|--------------|
| | Frialit® Internal Implant | 3.4 | M1.6 | YES | FRS 16 |
| | Frialit Internal Implant | 3.8 | M1.6 | YES | FRS 16 |
| | Frialit Internal Implant | 4.5 | M1.6 | YES | FRS 16 |
| | Frialit Internal Implant | 5.5 | M1.6 | YES | FRS 16 |
| | Frialit Internal Implant | 6.5 | M1.6 | YES | FRS 16 |
| CAMLOG® | CAMLOG Root Line | 3.8 | M1.6 | YES | FRS 16 |
| | CAMLOG Root Line | 4.3 | M1.6 | YES | FRS 16 |
| | CAMLOG Root Line | 5.0 | M2.0 | YES | FRS 20 |
| | CAMLOG Root Line | 6.0 | M2.0 | YES | FRS 20 |
| Biohorizons® | Maestro® External implants | 3.5 | M2.0 | Verify | FRS 20 |
| | Maestro External implants | 4.0 | M2.0 | Verify | FRS 20 |
| | Maestro External implants | 5.0 | M2.0 | Verify | FRS 20 |
| | Maestro External implants | 6.0 | M2.0 | Verify | FRS 20 |
| | Biohorizons Internal | 3.5 | 1-72UNF | Verify | FRS 172 |
| | Biohorizons Internal | 4.0 | 1-72UNF | Verify | FRS 172 |
| | Biohorizons Internal | 5.7 | 1-72UNF | Verify | FRS 172 |
| SPI® Implant | SPI Internal Implant | 3.5 | M1.4 | Verify | FRS 14 |
| | SPI Internal Implant | 4.5 | M1.6 | YES | FRS 16 |
| | SPI Internal Implant | 5.0 | M1.6 | YES | FRS 16 |
| | SPI Internal Implant | 6.0 | M1.6 | YES | FRS 16 |

ICL-IRK - 3/14-REV-3)

Distributed by:



4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 1-800-342-5454
 Outside The U.S.: +1-561-776-6700
 Fax: +1-561-776-1272
 www.biomet3i.com

Manufactured by:



e-space # 103, 104-1, 104-2, 105,
 106, 205, 212, 312, 509, 510, 511,
 10F, 36, 27 Gil, Digital-ro, Guro-gu,
 SEOUL, 152-789
 Rep. of Korea
 Tel.: 82-2-582-2885
 Fax: 82-2-582-2883
 www.neobiotech.co.kr
 Made in Korea

Branemark System, NobelActive, NobelReplace and NobelSpeedy are registered trademarks of Nobel Biocare. Camlog is a registered trademark of Camlog Biotechnologies AG. Nobel Biocare is a registered trademark of Nobel Biocare Services AG. Straumann is a registered trademark of Straumann Holding AG. Astra Tech, Frialit, Xive and Osseospeed are trademarks of Dentsply. NeoBiotech is a trademark of NEOBIOTECH CO., LTD. Advent, Screw-Vent, SwissPlus, Spine and Zimmer are registered trademarks of Zimmer Dental, Inc. Maestro and BioHorizons are registered trademarks of BioHorizons Implant Systems Inc. PrimaConnex, Renova and Stage-1 are registered trademarks of Keystone Dental, Inc. MIS and Seven are registered trademarks of Mis-Implants Technologies Ltd. SPI is a registered trademark of Thommen Medical. Certain is a registered trademark and Providing Solutions - One Patient At A Time and design are trademarks of BIOMET 3i LLC. BIOMET 3i and design are trademarks of BIOMET, Inc.
 © 2014 BIOMET 3i LLC. All rights reserved.