<b>DE</b>	Gebrauchsanweisung.....	2
<b>GB</b>	Instructions for use.....	14
<b>FR</b>	Mode d'emploi.....	26
<b>ES</b>	Instrucciones de uso.....	38
<b>IT</b>	Istruzioni per l'uso.....	50
<b>PT</b>	Instruções de uso.....	62
<b>NL</b>	Gebruiksaanwijzing.....	74
<b>SE</b>	Bruksanvisning.....	86
<b>DK</b>	Brugervejledning.....	98
<b>NO</b>	Bruksanvisning.....	110
<b>FI</b>	Käyttöohjeet.....	122

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Stand: 2023-03

**Medical Device – for use by dental health care professionals only.**

### **Intended purpose**

This dental accessory is suitable for patients requiring dental treatment for the following indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically, or they are accessories in the rehabilitation process.

RetraXil is an astringent retraction paste for the temporary displacement of the marginal gingiva and drying of the gingival sulcus. The paste, which contains aluminium chloride, is inserted into the gingival sulcus. In addition to the displacement effect of the paste, the astringent aluminium chloride ensures a haemostatic effect. Depending on the clinical situation and the practitioner's individual preference, the paste can be used as an alternative or in combination with retraction cords.

## **Indications**

For the temporary displacement of the marginal gingiva and drying of the gingival sulcus, e.g. for

- conventional or digital impressions,
- cementation of temporary and permanent restorations
- and the creation of Class II and V fillings.

## **Contraindications**

The use of this product is contraindicated in case of known or suspected allergies against components of this product. Do not use on patients with periodontal disease, furcation exposure or exposed bone.

## **Potential side effects**

Hypersensitivities to the product or its components cannot be excluded in individual cases. Local tissue damage if the application time is too long.

## **Product features**


### **Accessories to the packaging**

- A. Twisting aid (can be used multiple times after sterilisation)
- B. Application cannulas (for single use only – please discard after use)

### **Composition**

Phyllosilicates, water, aluminium chloride hexahydrate, mineral-based oleogel, vinyl siloxane, Pigment Blue 15:3

## Product application

### Preparation

- Thoroughly rinse the sulcus with water and gently air dry before using the paste. Remove and save the cap (Fig. 1).
- Depending on the clinical situation, a retraction cord can be inserted into the sulcus before the application of the retraction paste.
- Screw on one of the application cannulas supplied (Fig. 2), tighten it using the twisting aid provided and check that it is securely attached (Fig. 3).
- Protect the syringe with application cannula against contamination with a protective sleeve (Fig. 4).
- Apply a small amount of the retraction paste to a mixing pad and discard.

**Note:** Do not bend the cannula tip.

### Application


- Insert the application cannula into the sulcus and apply the retraction paste circularly around the tooth so that excess material wells out of the sulcus (Fig. 5).

**Note:** Do not use on prepared tooth surfaces – RetraXil can impair the adhesive bond of dentine-bonding agents.

- Leave the retraction paste in the sulcus for two minutes (Fig. 6).
  - No moisture should be allowed in during this time.

**Note:** A longer application time can lead to tissue damage.

### Rinse / remove

- Once the application time has ended, remove the retraction paste from the sulcus using an air-water mixture and a suction device; take care not to leave any residue (Fig. 7).
- Check the retraction result and remove any residue from the patient's sulcus/mouth.

**Note:** If RetraXil is not rinsed out of the mouth sufficiently, it can negatively affect the quality of the impression.

### Closure

- After application, remove the application cannula along with the protective sleeve using the twisting aid provided (Fig. 8)
- Place the syringe on a hygienically clean surface (Fig. 9)

- Discard the used application cannula along with the protective sleeve (Fig. 10) and reclose the syringe with the supplied cap using clean gloves (Fig. 11)

**Note:** When using syringes in mouth, use common hygienic protective sleeves to prevent contamination. Do not reuse after contamination.

- The twisting aid supplied must be sterilised after use (Fig. 12). For more information, see **“Cleaning, disinfection and sterilisation”**.

**Note:** Close RetraXil carefully after each use.

### **General information**

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

Visually inspect for damage before use. Damaged products must not be used. Do not reuse the syringes if there is evidence that the paste has been contaminated with body fluids.

### **Warnings and safety instructions**

- The product may irritate eyes. Avoid contact with eyes. Immediately and thoroughly rinse with plenty of water in case of contact with eyes and contact an ophthalmologist.
- The product may irritate skin. Avoid contact with skin. If the product comes into contact with the skin, thoroughly wash with soap and water immediately.
- Wear appropriate protective clothing, eye protection and gloves. The patient should use eye protection.
- The retraction paste must not be swallowed by the patient; this may result in an intolerance.
- Immediately clean paste residues from metal surfaces that have come into contact with the retraction paste, as they can corrode.

### **Storage conditions**

Use within ten weeks of opening the product. Do not store above 25°C (77°F). Store in a dry place. Do not use the material after the expiry date.

### **Disposal information**

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system.

European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

### **Cleaning and disinfection**

**WARNING: Please observe the country-specific regulations on hygiene standards and on the cleaning/disinfection of medical devices in dental practices. Wear suitable gloves. The twisting aids are supplied unsterilised. The twisting aid must be sterilised prior to first use and cleaned, disinfected and sterilised after each use.**

### **Processing restrictions**

Oxidising disinfectant solutions are not suitable. Strong acids and strong bases can damage the plastic.

## **TWISTING AID instructions**

### **In treatment room**

Unscrew the twisting aid from the syringe and discard the cannula along with the protective cover.

### **Retention and transport**

For transport outside of treatment room use a closed container. It is recommended that the twisting aid be reprocessed as soon as possible, but no more than two hours after use.

### **Preparing for decontamination**

Wear personal protective equipment (sturdy gloves, water-repellent protective gown, face mask or safety goggles and mask).

### **A) Cleaning and disinfection: by hand**

**Equipment:** non-protein-binding VAH-listed, virucidal disinfectant with cleaning action based on quaternary compounds (e.g. ID 212, Dürr Dental), lint-free cloths (e.g. Braun Wipes Eco), plastic brush (e.g. Interlock), if possible fully demineralised water (fully demineralised water, according to KRINKO/BfArM-recommendation free from facultative pathogenic microorganisms), tub for cleaning agents and disinfectants.

- 1) Prepare the cleaning solution according to the manufacturer's instructions (e.g. ID 212, 2%).
- 2) Immerse the twisting aid completely in the cleaning solution.
- 3) While immersed, brush the hard-to-reach areas of the immersed twisting aid with a soft brush. Pay attention to the critical hard-to-reach areas and cavities where it is not possible to visually assess the cleaning effect.
- 4) Immerse in the cleaning solution according to the manufacturer's instructions.

- 5) Remove the twisting aid from the cleaning solution and rinse thoroughly (e.g. 30 seconds) under running water.
- 6) Check for cleanliness; if contamination is still visible, repeat the above steps.
- 7) Prepare the disinfectant according to the manufacturer's instructions (e.g. ID 212, 2%).
- 8) Immerse the twisting aid completely in the disinfectant.
- 9) Immerse in the disinfectant according to the manufacturer's instructions (e.g. 5 minutes).
- 10) Remove the twisting aid from the disinfectant and rinse each one thoroughly (e.g. 30 seconds) with demineralised water.
- 11) Dry with a lint-free disposable cloth and/or with medical compressed air according to the European Pharmacopoeia.

Only validated processes may be used for the cleaning/disinfection of twisting aids.

## **B) Cleaning and disinfection: Automated WD [washer disinfectors]**

**Equipment:** cleaning and disinfection device according to DIN EN ISO 15883-1 and 2 with thermal program (temperature 90°C/194°F to 95°C/203°F),

**Cleaning agent:** mildly alkaline cleaning agent (e.g. neodisher® MediClean forte, Dr. Weigert).

- 1) Brush the hard-to-reach areas of the twisting aid for one minute under running potable water with a soft brush (e.g. Interlock). Pay attention to the critical hard-to-reach areas and cavities where it is not possible to assess the cleaning effect visually.
- 2) Place all twisting aids in a suitable mesh tray or place them on the load carrier so that all internal and external surfaces of the instruments are cleaned and disinfected.

<b>Program step</b>	<b>Water</b>	<b>Amounts required</b>	<b>Time</b>	<b>Temperature</b>
Pre-wash	CW (cold water)		5 min	
Cleaning solution dosing		Acc. to the manufacturer's specifications		Acc. to the manufacturer's specifications
Cleaning	DW (fully demineralised water)		10 min	55°C (131°F)
Rinsing	DW (fully demineralised water)		2 min	
Disinfection	DW (fully demineralised water)		3 min	A0 value > 3000 <sup>1</sup> (e.g. 90°C / 194°F, 5 min)
Drying			15 min	Up to 120°C (248°F)

<sup>1</sup> The authorities may issue other parameters for disinfection performance in their area of responsibility.

- 3) Remove all twisting aids at the end of the program.
- 4) Check that the load is dry and, if necessary, dry with medical compressed air in accordance with the European Pharmacopoeia or with a lint-free cloth.

- 5) After removal from the washer disinfecter, perform a visual check for cleanliness. If contamination is still visible, clean the medical devices by hand. Automated processing of the cleaned twisting aids must then follow.

### **Maintenance, inspection and testing**

All instruments are to be checked visually for cleanliness, integrity and functionality, if necessary, using an illuminated magnifying glass (3-6 dioptré).

All twisting aids must be checked for damage and wear. Damaged products must not be reused and must be disposed of.

## **C) Sterilisation**

**Device:** Steriliser according to DIN EN 285 or small steam steriliser according to DIN EN 13060, Type B process.

### **Packaging:**

- Equipment: Foil-paper packaging (e.g. steriCLIN, Item no. 3FKFS230110 and 3FKFS230112), sealing device (e.g. HAWO, Type 880 DC-V)  
A suitable procedure (sterile barrier system) must be used for packaging the twisting aids. Packaging according to DIN EN ISO 11607 (individual) or DIN 58953-9 (sets).
- Individual: A sterile barrier system (e.g. foil-paper packaging) according to DIN EN ISO 11607 should be used, one that is intended by the manufacturer for steam sterilisation. The packaging must be large enough so that the sealing seam is not under tension. (During the validation, double packaging is carried out in foil-paper packaging).

- Note: After the heat sealing process, the sealing seam must be checked visually for possible defects. In the event of errors, the packaging must be opened and the instrument repackaged and sealed.
- Sets: Sort twisting aids into the appropriate sieves so that they cannot be mechanically damaged during transport. Sterile barrier systems in accordance with EN ISO 11607 (e.g. sterilisation fleece) (the instruments can be single packed or double packed) or sterile goods containers in accordance with DIN 58953-9 are suitable for packaging the sieves. The weight of the sterile goods container must not exceed 10 kg.

**Procedure:** Steam sterilisation with fractional pre-vacuum, 134°C (273°F), holding time at least 3 min or 132°C (270°F) at least 4 min. Longer holding times are possible. In Germany, a holding time based on the KRINKO/BfArM-recommendation of at least 5 minutes is recommended.

- 1) Place the packaged twisting aid in the sterilisation chamber.
- 2) Start the program.
- 3) Remove the products at the end of the program and allow them to cool.

Check the packaging for any damage or moisture penetration. Compromised packaging is to be assessed as non-sterile. The instruments must be repackaged and sterilised.

**Storage:** Storage (protect from dust and moisture) and storage duration at the user's premises according to the specifications.

Only validated processes may be used for the sterilisation of twisting aids.

**Complaint notice**

Please quote batch number and article number or UDI in all correspondence about the product.

Article number **REF** and batch number **LOT** or UDI: See note on syringe/packaging.

According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

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Kulzer GmbH
Kulzer, LLC