INSTRUCTIONS FOR EC. STEAM 98 S

Single use pack for daily test of Bowie & Dick

EC. STEAM 98 S, thanks to its particular dimension, allows the operator to carry out a safe control of autoclave daily. Unlike the traditional devices, the use of the ready pack, presents significant advantages that guarantee to maintain the standards of proof for each test performed, associated to the elimination of complex practices of preparation of the same.

The ready single use pack **EC. STEAM 98 S** also presents important features for the safe reliability of the test:

- The packaging with big dimensions allows to collect immediately and precisely the presence of potential air bubbles inside the autoclave;
- The answer received remains unchanged in time;
- The external packaging has a changeable indicator that indicates at once if the pack was used or not;
- The material used for production are completely recyclable;
- Product in conformity with current regulations.

Product Description

Test Pack Bowie Dick **EC. STEAM 98 S** consists of: external packaging in printed carton with chemical indicator for steam cycles. The inside pack consists of 2 units of filtering barriers in cellulose. In the middle it is positioned the chemical indicator for the Bowie & Dick test.

The indicator sheet is printed with a chemical indicator sensitive to Temperature, Steam, Exposure Time to the previous variables that turn colour from yellow to dark grey/black.

Pack dimension 11,5x 11,5 in conformity to European regulations.

All materials used for the pack is suitable for medical use.

INSTRUCTIONS:

Before use: it is advisable to heat the autoclave carrying out an unladen sterilizing cycle.

After use: internal sheet in Bowie & Dick is not uniformly changed in grey/black highlighting a light bubble.

Result: the test has not been displayed to a right steam sterilization cycle, because the following minimum requirements:

- Exposure time: 3,5 minutes to a temperature of 134° C
- Loss or air infiltration caused by little function of chamber lining due to wear or breaking of the same.
- Loss in the sensor hatch connected to the chamber
- Pump for empty defective
- Cooling temperature of pump for empty too high
- Pressure sensor/switch defective
- Problems of condensation etc.

Suggestions: repeat the test. If the result is the same, suspend the use of the autoclave and ask for immediate after-sales service for a complete check. Fill the test report in every part and file it in the dedicated folder attaching the external sheet of Bowie & Dick packaging.

NO.1:

Before use: it is advisable to heat the autoclave carrying out an unladen sterilizing cycle.

After use: internal sheet of Bowie & Dick is uniformly changed in grey / light brown.

Result: Test has been exposed in a wrong way at an excessive temperature or to a lack of steam in the autoclave. Wet or overheated steam does not necessarily a failure of Bowie & Dick test, but can't be highlighted by the indicator sheet.

Suggestions: repeat the test. Fill the test report in every part and file it in the dedicated folder attaching the external sheet of Bowie & Dick packaging with process indicator.

NO. 2:

Before use: it is advisable to heat the autoclave carrying out an unladen sterilizing cycle.

After use: internal sheet of Bowie & Dick is uniformly changed in grey/black.

Result: Test has been exposed to saturated steam uniformly and correct and has been performed with success. It means that all base parameters, 134° C for 3,5 minutes, has been respected and the autoclave works under optimal conditions.

Suggestions: Fill the test report in every part and file it in the dedicated folder attaching the external sheet of Bowie & Dick packaging with the changed process indicator.

NOTE: a change of the internal sheet of Bowie & Dick from grey / intense black with burn signs indicate a complete and correct sterilization cycles, but with overexposure of all parameters: time, temperature, saturated steam. The test has been performed with success, but it is advisable to contact the after-sales service for a check of the autoclave, because an overexposed cycle can damage the autoclave and the material contained.

STORAGE OF DATA

Inside **EC. STEAM 98 S** pack, on the back of Bowie & Dick sheet, there is a complete card, easy to fill and practical to preserve, for storage of data.

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ECS Srl, Via Como, 71 - 23883 Brivio (LC) – Tel. 0341/681602 - Fax 0341/605485