Alminax Rite - Bite Index

Reinforced metallised wax shapes for bite registration

For professional use only.
Store in a cool, dry place away from sunlight

Manufactured by:
Associated Dental Products Limited
Purton, Swindon, Wiltshire, SN5 4HT, UK

Tel: (+44) 01793 770 256
Fax: (+44) 01793 772 256
Email: sales@kemdent.co.uk
Web: www.kemdent.co.uk


Fig. 1

Fig. 2

F009 Issue 1 01/10/10

Introduction

Alminax Rite-Bite Index gives a correct and reliable registration of hard tissue (tooth) contact ONLY and avoids soft tissue interference.

Both dentists and dental technicians find it easy to work with.

Technicians have noticed fewer problems with adjustments and corrections when their dentists use the Rite-Bite System.

Correct centric registration is the prerequisite for all prosthodontic work. Rite-Bite meets this criterion quickly and accurately.

Dr. Gunnar Edstrom: "Bite registration is one of the most important steps in crown and bridge work. Research has found more than 50% of all bite registrations coming to laboratories are wrong and unusable. A correct and reliable index is one of the most time saving factors in crown and bridge work".

Investigations have shown Kemdent Alminax Wax is the superior wax for bite registrations. It absorbs heat and softens quickly, yet cools easily and is very stable. Alminax Wax is the most suitable material for Rite-Bite Index.

Materials generally used for occlusal registrations often need to be stiffened. If whole sheets of wax, or wax supported by metal plates are used, the patient's tongue will be disturbed and this is likely to result in an inaccurate registration. The Rite-Bite Index has an integral reinforcement and a special cut-away shape so it does not interfere with the tongue - allowing the patient to bite comfortably in the correct position.

Instructions

1. Soften a Rite-Bite Index in water bath at 55°C (131°F).
2. Where there are preparations also soften a section of an unsupported index then press the two pieces together and melt the edges together with a hot wax knife. The different thickness of the wax does not matter because, heated in water, the wax will have a uniform softness.
3. Adjust the width of the index to the width of the patient's jaw (Fig. 1). The outer edge of the index should follow the buccal line of the teeth on the patient's largest jaw (usually the upper) Fig. 2. Never adjust a cool index.
4. Trim the length of the index to suit the patient's mouth, using scissors.
5. Re-soften the index in the water bath and place in the mouth so the outer edge follows the buccal line of the teeth on the larger jaw. It should be a little inside the outer cusp line. Make the registration. This uniform softness is the most important part of the bite registration because it allows the patient to bite easily and evenly into the Rite-Bite wax. If wax is heated over a flame it does not soften uniformly and this will cause the registration to be inaccurate. It is like trying to eat toffee by chewing straight up and down - the lower jaw automatically moves to the working side for the harder part.
6. Overbite: Use your finger to press the wax in the anterior region against the upper jaw.
7. Cool the index a little, but not completely, in the mouth, then take it out, turn it over to the smaller jaw side and cut out along the buccal cusps with scissors.
8. Cut away any wax distal to the last molar and also wax which has pressed against the soft tissue because the hard model will "hang" on these areas and give an inaccurate position.
9. Warm the index a little in warm water, if necessary, replace in the mouth and check the registration is perfect.
10. Cool the index completely, remove it from the mouth.
11. Disinfect the wax and pack it in a padded box for safe delivery to the laboratory.

Visit www.kemdent.com for our full range of
Cross - Infection Control Products:
PracticeSafe, ChairSafe & InstrumentSafe Disinfectants

Product Information

Alminax Rite-Bite Index ready-shaped to mouth arch and reinforced round inner edge. Supplied in a box of 30, containing 24 reinforced and 6 plain arches.

Product Code: DWS211

- 2.3mm thick
- Full arch bite wafers
- Preformed to U shape
- Colour: Silver Grey
- Softens at 55°C
- Reinforced with metal strip

Full range of Alminax products available:

Alminax Bite Wax sheets - 143mm x 73mm x 2.6mm thick.

Product Codes: DWS219 - 250g box, DWS220 - 500g box.

Alminax Bite Wax Full-arch - ready shaped to mouth arch.

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