

Product information

SENDOLINE AB

K- File

The K-Files is designed with four cutting edges and a 60° tip with 135° cutting angle with good cutting ability. This together with its tight spiral ensure a safe preparation.

Use the K-File with a twisting/filing movement.

The K-File is also available in nickel titanium named [NiTi K-File](#).

ISO/Sizes: 06-140

Length: 21 mm, 25 mm, 28 mm, 31 mm

Package: Individual or assorted (in sizes 15-40 or 45-80) in blister. All packages come with 6 pcs/blister

Material: Operating part: Stainless steel. Handle: Plastic. Stopper: Silicone.

