


# Fit Strip™

**F**lexible  
**I**nterproximal finishers with a  
**T**wist



**FitStrips** are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison Dental Solutions.**

**One handy handle!** FitStrips' simple, detachable handle makes a world of difference for both clinician ease of use and patient comfort. Your hand stays outside the mouth, granting a clear field of view and easier operation.

**Do the Twist!** Twist the color-coded barrel to adjust FitStrips' curvature and it's automatically locked in place. It doesn't get any easier than that!



**FitStrips make your procedures easier!**

- Composite restorations: contouring/finishing all interproximal surfaces
- Orthodontics: IPR/tooth slenderizing
- Crown and bridge: cement removal and clean-up

**Ready for a marathon!** FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.

# Fit Strip™

**F**lexible  
**I**nterproximal finishers with a  
**T**wist



### FPSK01 Kit contains:

- 2 handles
- 2 serrated FitStrips
- 1 single-sided FitStrip of each grit (4 pcs)
- 1 double-sided FitStrip of each grit (4 pcs)

### Variety is the key to happiness!

One FitStrips starter kit contains everything you need:

- Single-sided diamond strips
- Double-sided diamond strips
- Multiple diamond grits
- Even serrated FitStrips are included!






## ADDITIONAL KITS

All refills come in a package of 10 FitStrips with 2 handles.

Single Sided Kit 2 ea of all 4 grits plus 2 serrated	FPSK02
Double Side Kit 2 ea of all 4 grits plus 2 serrated	FPSK03

## REFILLS

All refills come in a package of 10 FitStrips with 2 handles.

		Single Sided	Double Sided
<b>Super fine</b> 15 Micron 0.11 Thick		FPSXFSS	FPSXFDS
<b>Fine</b> 30 Micron 0.15 Thick		FPSFSS	FPSFDS
<b>Medium</b> 46 Micron 0.21 Thick		FPSMSS	FPSMDS
<b>Coarse</b> 90 Micron 0.30 Thick		FPSCSS	FPSCDS
<b>Serrated</b> 0.05 Thick		FPSSAW	

For more information go to [www.garrisdental.com](http://www.garrisdental.com). Be sure to check out the website to see the full line of innovative Garrison products and our 6-month money-back guarantee.

**Garrison**  
Dental Solutions

150 DeWitt Lane  
Spring Lake, MI USA  
49456

To order: Call toll free **888.437.0032**  
or online at [www.garrisdental.com](http://www.garrisdental.com)



Mention Code  
**IN1016FS**