

# Granulation Tissue Instruments

acc. to Dr. Mike Sabeti, DDS  
University of California, San Francisco



**3305**

2 mm diameter,  
short and long  
PEEK handle, yellow



**3306**

4 mm diameter,  
short and long  
PEEK handle, green

This surgical device is capable of efficiently and effectively removing the granulation tissue encountered in periodontal surgery, endodontic surgery and extraction socket.

The current surgical devices currently utilized in the surgical removal of granulation tissue are extremely time consuming. The granulation tissue devices have been designed to effectively remove granulation tissue. The proposed surgical device consists of a loop in the middle with sharp edges that grab and cut the granulation tissue fast and more efficient.

