



PRESIDENT[®]
Impression Tray System


**If they weren't transparent you'd think
they were made of steel!**

PRESIDENT®

Impression Tray System






Rigid by design



Anatomically correct



Built-in retention



Feature	Benefit
Rigid by design	A patented system of intersecting ribs and rails provides unprecedented rigidity in a plastic tray
Anatomically correct	A statistical database accumulated over 20 years has enabled Coltène / Whaledent to define the ideal form for each tray system
Built-in retention	The President Tray System's design aggressively retains all types of impression material.
Transparent	Easy and precise positioning of the tray.
Modifiable form	The tray can individually be adjusted with heat.
Single use	No risk of contamination, no time-consuming cleaning.
Economical	Good value for money

Product Ordering Information: 1 Bag of 12 Trays

Description	Dimensions		
6901 1 Large Upper	A 83.1 mm/3.270"	B 63.2 mm/2.490"	C 18.7 mm/.738"
6902 2 Large Lower	A 79.9 mm/3.147"	B 62.2 mm/2.450"	C 20.0 mm/.788"
6903 3 Medium Upper	A 75.4 mm/2.968"	B 58.3 mm/2.295"	C 17.3 mm/.683"
6904 4 Medium Lower	A 69.5 mm/2.738"	B 58.5 mm/2.303"	C 19.3 mm/.759"
6905 5 Small Upper	A 69.6 mm/2.714"	B 51.7 mm/2.036"	C 17.1 mm/.671"
6906 6 Small Lower	A 61.7 mm/2.431"	B 50.5 mm/1.990"	C 18.3 mm/.720"
6907 7 Quadrant UL/LR	A 47.1 mm/1.854"	B 58.5 mm/2.303"	C 19.3 mm/.759"
6908 8 Quadrant UR/LL	A 47.1 mm/1.854"	B 58.5 mm/2.303"	C 19.3 mm/.759"
6909 9 Anterior	A 65.3 mm/2.569"	B 39.4 mm/1.554"	C 19.2 mm/.756"

© 2013 Coltène/Whaledent AG – www.coltene.com

Coltène/Whaledent AG

Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel +41 71 757 5300
Fax +41 71 757 5301
info.ch@coltene.com

Coltène/Whaledent Inc.

235 Ascot Parkway
Cuyahoga Falls, Ohio 44223 / USA
Tel +1 330 916 8800
Fax +1 330 916 7077
info.us@coltene.com

Coltène/Whaledent Ltd.

The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF / UK
Tel +44 1444 235486
Fax +44 1444 870640
info.uk@coltene.com