



# Temporary success

**NOBETEC<sup>®</sup>**  
ZOE-cement

“When nothing else works!”

NOBETEC®



# NOBETEC® – ZOE-cement

NOBETEC® is a tissue compatible multi-purpose ZOE-cement. The powder and liquid offers the possibility to mix different consistencies for various applications.

## COMPOSITION

Powder:  
Zinc oxide 60%  
Colofonium 40%

Liquid:  
Eugenol 60%,  
Balsam of Canada 35%  
Balsam of Peru 5%

NOBETEC® is known for its strength!

## CEMENTATION OF TEMPORARY CROWNS AND BRIDGES

It is also easy to remove without falling off ahead of schedule when used for cementation of temporary crowns and bridges.

## TEMPORARY FILLING

When using NOBETEC as temporary filling the adhesion is excellent but at the same time it is easily removed by means of an excavator or a burr drill.

## SURGICAL DRESSING

NOBETEC protects wound surfaces, is tissue compatible and encourages rapid epithelisation which normally occurs within 2 weeks.

| ● FACT                             | ● BENEFIT                            |
|------------------------------------|--------------------------------------|
| Multi purpose                      | Three indications – only one product |
| Eugenol, Balsam of Canada and Peru | Pain relieving and soothing capacity |
| Excellent adhesion                 | Stays in place as long as needed     |

## DID YOU KNOW?

Thanks to the total etch technique any remaining eugenol is removed and does not affect the bonding capacity. Eugenol, Balsam of Peru and Balsam of Canada have all pain relieving and soothing capacity.

## PACKAGING:

Powder: 100 gram in a plastic container  
Liquid: 20 ml in a glass bottle



## Different personalities – One happy family



Nordiska Dental develops, manufactures and markets dental materials with packs and accessories that make dentistry simpler and safer.

