

BEAUTIFUL RETRACTION

The image shows a white box for the RetraXil Starter Pack and a syringe. The box features the RetraXil logo, the text 'Starter Pack', and two lines of German text: 'Adstringierende Retraktionspaste' and 'Adstringent retraction paste'. The syringe is blue and white with the RetraXil logo on it. The Kulzer logo is also visible on the box.

RetraXil®

Starter Pack
Adstringierende Retraktionspaste
Adstringent retraction paste

 KULZER
MITSUI CHEMICALS GROUP

RetraXil®

Good sulcus exposure – minimally invasive application

RetraXil® provides a cleanly opened gingival sulcus and ensures beautifully exposed, blood-free preparation margins. Its syringe sits comfortably in the hand and supports precise dosing of the material.

Your benefits:

- » Atraumatic application
- » Adequate tissue management
- » Multidose stable syringe
- » No curls, no tear-off, non-stick

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

RetraXil® – Beautiful Retraction

Efficient, gingival retraction with highest standards in handling and safety

Safety during application

- No risk of systemic overdose
- Clear colour contrast to the gingiva and teeth offers high visibility and good control when rinsing
- Twisting aid offers secure fit of the cannula on the syringe and enables the cannula to be twisted on and off safely

Good handling features

- Thin cannula offers precise and atraumatic application into the sulcus and interproximal areas
- Compatible with analogue and digital impression techniques

Economical and time-saving

- Fewer work steps compared to retraction cord technique



The consistency of a paste – the stability of a cord

RetraXil® offers atraumatic and easy access to the sulcus thanks to the thinnest cannula on the market

RetraXil®	Competitor 1	Competitor 2	Competitor 3	Competitor 4	Competitor 5
Ø 0.9 mm	Ø 1.25 mm	Ø 1.25 mm	Ø 1.4 mm	Ø 1.4 mm	Ø 1.6 mm
Cannula	PLT	Cannula	PLT	PLT	Cannula

Art. Code	Article	Content
Products		
66098180	RetraXil Starter Pack	4 x 1 g RetraXil 2 x Twisting aid 20 x Cannula
66098179	RetraXil Value Pack	20 x 1 g RetraXil 10 x Twisting aid 100 x Cannula
Refills		
66087214	RetraXil Twisting Aid Refill	2 x 10 Twisting aids
66087215	RetraXil Cannula Refill	5 x 20 Cannulae



kulzer.com/retraxil