

As easy as 1-2-3



Repair Kit from GC

Smart products
and guidelines
for ALL your intra-oral
repair procedures



GC

The GC Repair Kit, chairside repair made easy...

A **complete set** of all products and guidelines needed to cover all your **routine intra-oral repairs** in an **easy and smart way**.



Metal-based indirect restorations
(precious and non-precious)



Indirect restorations made of
glass-ceramics, zirconia, alumina or
hybrid ceramics (e.g. CERASMART)



Composite restorations
(direct or indirect)

What's in it for you? Three components to fit all repairs!

Ceramic Primer II

Following the zero compromise approach, we suggest the use of a **separate silane primer** (Ceramic Primer II) together with our adhesive G-Premio BOND, for a durable adhesion to **glass ceramics** and **hybrid ceramics**.



GRADIA PLUS Opaque

This **light-cured opaquer** efficiently masks dark backgrounds, to ensure high aesthetics in case of porcelain fused to metal restorations.

G-Premio BOND

A **universal bonding** which offers **optimal adhesion to all substrates** thanks to 3 different monomers.

4-MET

for a reliable
adhesion to **enamel**
& **dentin**

MDP

for an excellent
adhesion to **enamel**
& **dentin, zirconia,**
alumina and
non-precious
metals

MDTP

to bond to
precious metals*



For the step-by-step, consult
also the GC Restorative
Dentistry Guides!

Article Number	Description
901139	GC Repair Kit: Ceramic Primer II (5 unitdoses), G-Premio BOND (5 unitdoses), GRADIA PLUS OAI (1 syringe 2.0mL)



*MDTP:
Methacryloyloxydecyl dihydrogen
thiophosphate to bond to precious metals

GC EUROPE N.V.

Head Office
Researchpark
Haasrode-Leuven 1240
Interleuvenlaan 33
B-3001 Leuven
Tel. +32.16.74.10.00
Fax. +32.16.40.48.32
info@gceurope.com
http://www.gceurope.com

GC UNITED KINGDOM Ltd.

16-23, Coopers Court
Newport Pagnell
UK-Bucks. MK16 8JS
Tel. +44.1908.218.999
Fax. +44.1908.218.900
info@uk.gceurope.com
http://uk.gceurope.com