



# QLICKSMART®

## Safety Catalogue 2025

# QLICKSMART®

... and be safe

## BladeFLASK



"Very convenient and safe to remove scalpel blades..."  
-Lab Technician

"The nursing and medical staff feel safer with the Qlicksmart BladeFLASK"  
-Dr. Daryl Wall

"Thanks to whoever invented the Qlicksmart BladeFLASK. Sheer genius."  
-Mark Kosinski

"We feel much safer using them. Thank you so much!"  
-Angela Connolly, MLT

<input checked="" type="checkbox"/>	Single-handed blade removal	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	Removes and contains 100 blades	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	Sharps container per AS 4031:1992	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	Can be mounted with Universal Mounting Bracket	<input checked="" type="checkbox"/>
<b>Compatible with:</b>		
<input type="checkbox"/>	All commercially available scalpel handles and blades	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	#3, 3L, 4, 4L, 7, 9, Lawton long handles	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	#10, 10a, 11, 12, 13 - 25a, 27, and 36 blades	<input checked="" type="checkbox"/>
<input type="checkbox"/>	Bulbous handles (PM range)	<input checked="" type="checkbox"/>
<input type="checkbox"/>	Hexagonal handles	<input checked="" type="checkbox"/>
<input type="checkbox"/>	Circular handles	<input checked="" type="checkbox"/>
<input type="checkbox"/>	Barron handles	<input checked="" type="checkbox"/>

## BladeFlask EVO



**Bracket**  
Enables single-handed blade removal

The Qlicksmart BladeFlask EVO can accommodate multiple types of #3 & #4 fitment handles. This includes large autopsy scalpels, hexagonal, Barron, bulbous, thick and round handles, and the new ergonomic handles

Use **SnapIT** and **CheckCLIP** together to protect against...



**Injury**



**Contamination**



**Error**



**Ampoule opener**



**Medication Identification**



**Safety**

# BladeFLASK™

## scalpel blade remover & sharps container

Remove, contain, and dispose of scalpel blades in one click with the BladeFLASK™. The world's pioneering single-handed scalpel blade remover enables healthcare professionals to quickly and safely remove up to 100 scalpel blades.

- ✓ Puncture-proof, compliant with UN recommendations for handling dangerous waste
- ✓ Audible “click” to signify scalpel blade removal
- ✓ Built-in counter with automatic shutoff when maximum capacity is reached
- ✓ OSHA compliant with single-handed activation when used with the reusable Mounting Bracket
- ✓ Mount on a wall or desk with the reusable Mounting Bracket
- ✓ Compatible with most brands of general surgery blades and flat handles

Video Instructions

Scan QR Code



Extra safety when used with the Quicksmart Mounting Bracket, which affixes the BladeFLASK to the wall or bench to allow for single-handed scalpel blade removal



BladeFLASK Product Code	QS1000
Quicksmart Mounting Bracket Product Code (Sold separately)	QS2000
Compatible scalpel blades	#10, 10a, 11, 11P, 12, 14, 15, 15a, 15c, 16, 19, 20, 21, 22, 22a, 23, 24, 25, 25a, 27, 36
Compatible scalpel handles	#3, 3L, 4, 4L, 5, 7, 9, Lawton long handle
Incompatible scalpel handles	Baron handle, Beaver handle, disposable handles

# BladeFlask EVO™

## scalpel blade remover

The Qlicksmart BladeFlask EVO™ can accommodate even more types of scalpel handles and blades including large autopsy scalpels, thick and round handles, even hexagonal shaped versions - giving you greater freedom in using any scalpel you want. It is so simple to use with just one hand. Embrace safety with the new standard in scalpel blade removals.

- ✓ Puncture-proof, compliant with UN recommendations for handling dangerous waste
- ✓ OSHA compliant with single-handed activation when used with the reusable Mounting Bracket
- ✓ Built-in counter with automatic shutoff when maximum capacity of 100 blades is reached
- ✓ Compatible with all brands of general surgery blades including large autopsy scalpels, thick and round handles, even hexagonal shaped versions
- ✓ Mount on a wall or desk with the reusable Mounting Bracket
- ✓ Ideal for mortuaries, plastic surgeries autopsy stations and CSSD

Video Instructions

Scan QR Code



Extra safety when used with the Qlicksmart Mounting Bracket, which affixes the BladeFLASK to the wall or bench to allow for single-handed scalpel blade removal



BladeFLASK EVO Product Code	QS3000
Qlicksmart Mounting Bracket Product Code (Sold separately)	QS2000
Compatible scalpel blades	All commercially available BP scalpel blades including (but not limited to) up to size PM60A, PM60B, PM70.
Compatible scalpel handles	All commercially available scalpel handles including (but not limited to) PM6, PM8.

# Single Use Range

## sterile scalpel blade removers + needle counter

Quicksmart's range of single-use scalpel blade removers are engineered for the unique needs of the Operating Theatre. They utilise patented automatic safety technology for safe and efficient sharps removal, containment, counting, and disposal.

- ✓ Audible “click” to signify scalpel blade removal
- ✓ OSHA compliant with single-handed activation
- ✓ Transparent cartridges for easy containment and counting
- ✓ Available off-the-shelf (sterile) or as part of a surgical procedure pack
- ✓ Red arrows and guides show proper insert location
- ✓ BladeNeedleSYSTEM™ includes a 20-count foam needle counter with an Overlapping lid and magnetic surface



Video Instructions

Scan QR Code



Video Instructions

Scan QR Code



Video Instructions

Scan QR Code



	BladeCASSETTE	BladeNeedleSYSTEM	BladeSINGLE
Product Code	QSCS3100	QSBNSS	
Blade Capacity	3	2	1
Compatible scalpel blades	#10, 10a, 11, 11p, 12, 14, 15, 15a, 15c, 16, 19, 20, 21, 22, 22a, 23, 24, 25, 25a, 27 36, 6, 9, E11, Sabre D/15	#10, 10a, 11, 11p, 12, 14, 15, 15a, 15c, 16, 19, 20, 21, 22, 22a, 23, 24, 25, 25a, 27 36, 6, 9, E11, Sabre D/15	#10, 10a, 11, 11p, 12, 14, 15, 15a, 15c, 16, 19, 20, 21, 22, 22a, 23, 24, 25, 25a, 27 36, 6, 9, E11, Sabre D/15
Compatible scalpel handles	#3, 3L, 4, 4L, 5, 7, 9, Baron, Lawton long handle, cylindrical handle	#3, 3L, 4, 4L, 5, 7, 9, Baron, Lawton long handle, cylindrical handle	#3, 3L, 4, 4L, 5, 7, 9, Baron, Lawton long handle, cylindrical handle
Incompatible scalpel handles	Beaver handle, Disposable handles	Beaver handle, Disposable handles	Beaver handle, Disposable handles

# SnapIT™

## reusable ampoule openers

Quicksmart's reusable range of SnapIT™ ampoule openers allow healthcare professionals to break open glass ampoules swiftly and without anxiety. Three models of SnapIT™ have been designed to fit the needs of individuals or organisations.

- ✓ Easy “insert and snap” action
- ✓ Prevents injuries or medication spillage by facilitating a clean break
- ✓ SnapIT Personal and SnapIT Trolley editions can be autoclaved to 130° Celsius
- ✓ Can be disassembled for cleaning
- ✓ Keyring attached for safe ampoule opening on-the-go
- ✓ SnapIT Personal and SnapIT Lite editions most used ampoule openers

Video Instructions

Scan QR Code



	SnapIT Personal	SnapIT Lite		SnapIT Trolley		
Product Code	QSSN-PR	QSSN-LITE-R	QSSN-LITE-L	QSSN-TR	QSSN-TL	QSSN-TXL
Material	Aluminium	Lightweight Plastic		Aluminium		
Available sizes	Regular	Regular	Large	Regular	Large	Extra Large
Compatible ampoule size (approximate)	1-15mL	1-15mL	5-25mL	1-15mL	5-25mL	20-30mL
Colours	blue	grey	burgundy	black	red	blue

# CheckCLIP™

## medication identification

Qlicksmart's CheckCLIP™ medication identification clips protect patients by helping to prevent medication errors. Invented by paramedics, the CheckCLIP™ provide accurate and efficient medication identification during busy and stressful times.

- ✓ Attaches to Luer and Luer-lock syringes
- ✓ Doesn't obstruct dosage markings on the syringe
- ✓ Available for vials or ampoules
- ✓ Barbs on clip prevent detachment
- ✓ Allows both the medication and syringe to be held in your hands when drawing up
- ✓ Bright colours for high visibility

Video Instructions

Scan QR Code



	CheckCLIP for Ampoules	CheckCLIP for Vials
Product Code	QSCC100	QSCV100
Compatible Syringe Size	1-20mL Luer and Luer-Lock	1-20mL Luer and Luer-Lock
Compatible Medication Size (Approximate)	1-10mL glass and plastic ampoules	Vials with diameter up to 28mm
Colour	pink	green

**Wilt u snij- en prikincidenten,  
huidverwondingen en infecties voorkomen?**

Gebruik dan Qlicksmart producten.  
Veiligheid voor personeel en patiënten  
is onze prioriteit.

**100% veilig, hygiënisch, preventief  
en kostenbesparend!**



Qlicksmart voldoet aan de nieuwe Wet- en Regelgeving werken met Scherpe Medische Hulpmiddelen art.4.97 punt 2c. Nadere informatie en bestellingen bij uw reguliere groothandel of Medzorg Nederland B.V.


**Importer / Distributor**

Medzorg Nederland BV  
Satellietbaan 14 B  
2181 MH Hillegom



 0031 (0)252 536200

 [www.qlicksmart.nl](http://www.qlicksmart.nl)

 [orders@medzorg.nl](mailto:orders@medzorg.nl)  
[info@qlicksmart.nl](mailto:info@qlicksmart.nl)

**Sub-Distributor**

