




ART³: A NEW WAY TO BULK FILL

Achieve lifelike restorations with the next generation of clinically proven Adaptive Response Technology (ART) found in Harmonize™ and OptiShade Universal Composite. Powered by Enhanced Adaptive Response Technology, SimpliShade Bulk Fill offers easy handling and low shrinkage.




Adaptive Rheological Response Technology for Ease of Sculpting

ART³ uses rheological modifiers that reduce the composite's viscosity as you manipulate it. When you're done shaping, viscosity recovers and anatomy stabilizes before light-curing — reducing the need for resculpting.

Adaptive Color-Blending Response Technology for Beautiful Blending

Achieve exceptional blending with special nanofillers that mimic enamel structures and interact with light in a natural way. The restoration and tooth become virtually indistinguishable.




Human Enamel

SimpliShade Bulk Fill

Adaptive Shrinkage Stress Response Technology for Durable Results

During polymerization, SimpliShade Bulk Fill uses low shrinkage resin and innovative nanofiller design, resin-infused nanofillers boost elasticity for stress relief and a special coating eases interfacial stress. Uniting three shrinkage reduction mechanisms, this enhances marginal adaptation and integrity by minimizing composite restoration stress on the adhesive layer for enhanced marginal adaptation.




Learn more about ART³ and the effect on shrinkage stress

Ordering Information



SimpliShade Bulk Fill Flow

Part #	Item Description
37217	SimpliShade Bulk Fill Flow Syringe
37218	SimpliShade Bulk Fill Flow Unidose 20 Pack



KEEP IT SIMPLE

Choose SimpliShade Bulk Fill or SimpliShade™ Universal Composite for your capping layer.

1. SimpliShade Bulk Fill Flow 2. Capping layer with SimpliShade Bulk Fill or SimpliShade Universal Composite



SimpliShade Bulk Fill

Part #	Item Description
37211	SimpliShade Bulk Fill Composite Syringe
37212	SimpliShade Bulk Fill Composite Unidose 20 Pack



SimpliShade Universal Composite

Unidose		Syringe	
Part #	Item Description	Part #	Item Description
37006	Light - 20 Pack	37001	Light
37007	Medium - 20 Pack	37002	Medium
37008	Dark - 10 Pack	37003	Dark
37009	Bleach White - 10 Pack	37004	Bleach White
37010	Universal Opaque - 10 Pack	37005	Universal Opaque



For an ultra high shine finish, use Kerr's line of Polishers!



One Shade. One Step. A Better Way to Bulk Fill.

For more information scan the QR Code or visit go.kerrdental.com/simplishadebulkfill

SimpliShade™ BULK FILL FLOW



Single-Shade Simplicity Meets Superior Adaptability

Reach every area of the tooth's anatomy while saving time with SimpliShade Bulk Fill Flow. Requiring just one shade to match all 16 VITA® classical shades, SimpliShade Bulk Fill Flow easily adapts to the cavity floor — with exceptional durability and depth of cure.

kerrdental.com

©2023 Kerr Corporation. All trademarks applicable to Kerr and/or its products are property of Kerr Corporation. All other trademarks are property of their respective owners. All Rights Reserved. Data on File. MKT-22-0880_UK Rev 0

Kerr™

Kerr™

ONE SHADE IS ALL IT TAKES

Get the adaptability, durability, depth of cure, and simplicity you need with SimpliShade Bulk Fill Flow.

SimpliShade™ Bulk Fill Flow is indicated for base/liner under **Class I** and **Class II** direct restorations, pit and fissure sealant, restoration of small cavities, **Class III** and **Class V** restorations, blocking out undercuts, small core build-ups, and repair of resin and acrylic provisional restorations.

Flowable Features Single-Shade Simplicity



One-Shade System

Using one shade that blends seamlessly to the cavity walls with all 16 VITA® classical shades, SimpliShade Bulk Fill Flow saves you time selecting shades and managing inventory.



Save Time

Simplify your workflow with just one step to cure up to 4mm, before adding a capping layer.



Go with the Flow

Fill every area of the tooth's anatomy with a self-leveling bulk fill flowable that easily adapts to the cavity floor.



Lower Shrinkage

Helps assure a long-lasting restoration with lower polymerization shrinkage and shrinkage stress when compared to other bulk fill brands.




Add Confidence

Restore confidently with a flowable bulk fill composite that offers comparable or higher mechanical strength than the competitors.


Shrinkage Stress

SimpliShade Bulk Fill Flow offers lower shrinkage stress than Competitor 1, Competitor 3, Competitor 2, and Competitor 4, resulting in lower marginal degradation and postoperative sensitivity for a long-lasting restoration.




Mechanical Strength

SimpliShade Bulk Fill Flow has comparable or higher mechanical strength than Competitor 5, Competitor 1 and Competitor 2, resulting in better resistance to chipping and fracture over time.



Radiopacity (AI%)




Radiopacity

Tests show that SimpliShade Bulk Fill Flow provides higher radiopacity compared to Competitor 1, Competitor 3, Competitor 4. Higher radiopacity enables clinicians to easily identify restorations radiographically.

Wear

SimpliShade Bulk Fill Flow offers a lower wear rate at interproximal contacts than Competitor 5. Low wear rate results in a long-lasting restoration.



Adaptation

Easily adaptable, SimpliShade Bulk Fill Flow fills every area of the tooth's anatomy.

